

# Customize Your Case


Tailored Ankle Fusion Configurations

**EPIC**  
EXTREMITY





Designed for ankle arthrodesis, the EPIC Extremity Ankle Fusion Plate System provides surgeons with versatile treatment options in a customizable tray. Surgeons can choose to use either cannulated screws, plates or a combination of cannulated screws and plates based on their patient's specific needs.



## KEY FEATURES

- Complete anatomic titanium plate portfolio
- Right- and Left-Specific
- Multiple sizes and lengths
- Compression slots
- Low-profile plates with tapering thickness, max thickness at 3.5mm

## ANATOMIC FUSION PLATE OPTIONS

- A. Anterior: Standard and small tension band
- B. Anterolateral: Four options – Left, right, standard and long
- C. Lateral: Four options – TT, TTC standard and long
- D. Posterior: Two options – TT and TTC

## VARIABLE LOCKING AND NON-LOCKING SCREW OPTIONS

- 4.5mm variable locking and non-locking screws in lengths 16-40mm in 2mm increments, 45mm and 50mm
- 5.5mm variable locking and non-locking screws in lengths 26-50mm in 2mm increments, 55mm and 60mm

